

Surgicare®

DUAL PREMIUM DHD

**Sterile Latex Surgical Powderfree Gloves
Damp Hand Donning with Moisturizer**

Surgicare Dual Premium DHD Gloves are made from Natural Rubber Latex. These gloves are anatomically shaped with thumb ball effect, reducing hand fatigue during surgical procedures. Surgicare Dual provides primary and secondary protection for both doctor and patient during high risk surgeries. These gloves are non-chlorinated.

FEATURES

- ▶ **Damp Hand Donning (DHD) with moisturizer.** DHD Technology assists in donning gloves over damp hands and moisturizer keeps the skin hydrated & soft.
- ▶ **Visual warning of a puncture on secondary glove in the presence of fluid, as a green stain appears.**
- ▶ Recommended for High Risk, Orthopaedic & Trauma surgery.
- ▶ Protects against pathogen transmission and is a barrier to cross infection.
- ▶ **The powderfree Under Glove protects the user against latex protein allergies and cross contamination from glove tear/puncture.**
- ▶ **The green colour of the 'Underglove' helps the Surgeon quickly recognize outer glove perforation.**
- ▶ Inner Polymer Coating provides smooth and easy donning of gloves. Specially designed mould produces gloves of high comfort.
- ▶ Protein content no more than 50µg of protein per gram & dm².
- ▶ Manufactured in dust free environment with low bioburden and therefore reduced pyrogen.
- ▶ Latex Film structure passes Viral Penetration Test (ASTM F1671)

STANDARDS

- ▶ ASTM : D-3577 & EN 455-1, EN 455-2, EN 455-3 & ISO 10282
- ▶ CE Mark - MDR: Regulation (EU) 2017/745, Class IIa.
- ▶ Passes Biocompatibility test as per ISO 10993-10
- ▶ Outer Glove Tested against Chemotherapy drugs as per ASTM D6978.

Recommended for:

- ▶ **High Risk, Orthopaedic & Trauma Surgery**



Kanam
KANAM LATEX
INDUSTRIES PVT. LTD.

The trusted name the world over

www.kanamlatex.com



Ooppotttil Buildings, K.K. Road
Kottayam-686 001, Kerala, India
Tel: 0481-2560108, 2563513
Fax: 0481-2563614
E-mail: customercare@kanamlatex.com
ho@kanamlatex.com

AN ISO 9001:2015, ISO 13485:2016 & ISO 14001:2015 CERTIFIED COMPANY

This product is CE Certified