



Powdered Non Sterile Latex Examination Gloves

Kaltex Plus Powdered Non Sterile Latex Examination Gloves are made from Natural Rubber Latex. These gloves provide excellent barrier protection for both healthcare worker and patient.

SAFETY FEATURES

- ▶ Enhanced elasticity, tear strength and donning.
- ▶ Beaded cuff with minimum length of 230 mm. Can offer 290 mm length also.
- ▶ Powdered with Modified Bio-Absorbable Corn starch, to prevent wound infections and granulomas.

STANDARDS

- ▶ ASTM D 3578.
- ▶ Registered with USFDA (510K).

QUALITY CONTROL

- ▶ Factory standard of 1.5 AQL for pinholes.



Kanam
KANAM LATEX
INDUSTRIES PVT. LTD.

The trusted name the world over

Ooppottil Buildings, K.K. Road
Kottayam-686 001, Kerala, India
Tel: 0481-2560108, 2563513
Fax: 0481-2563614

E-mail: customercare@kanamlatex.com
ho@kanamlatex.com

AN ISO 9001:2015, ISO 13485:2016 & ISO 14001:2015 CERTIFIED COMPANY

www.kanamlatex.com

Product Information

TYPE

Powdered Non Sterile Latex Examination Gloves

PRIMARY MATERIAL

Natural Rubber Latex (NR).

COLOUR

Creamy White.

DESIGN

Ambidextrous, Straight Finger, Beaded Cuff, Smooth surface.

POWDER AMOUNT

<10mg/dm².

SHELF LIFE

3 years in original packaging if stored as per storage instructions.



PHYSICAL DIMENSIONS

Size	Minimum Length (mm)	Width (mm)	Thickness (mm) Min	
			Palm	Finger
XS	230	70 ± 10	0.08	0.08
S	230	80 ± 10	0.08	0.08
M	230	95 ± 10	0.08	0.08
L	230	111 ± 10	0.08	0.08
XL	230	120 ± 10	0.08	0.08

PHYSICAL PROPERTIES

Characteristics	Before Ageing	After Ageing 100±2°C for 22±0.3hrs
Tensile Strength (Mpa) min.	18	14
Ultimate Elongation (%) min.	650	500
Stress at 500% Elongation (Mpa)	5.5 max.	N/A

PRE-SHIPMENT QUALITY INSPECTION

Sampling Procedure: ISO 2859 Part I

Sampling Plan: Single sampling, Normal Inspection

Characteristics	Inspection Level	AQL
Freedom from Holes	G-1	1.5
Visual attributes – Major Defects	G-1	1.5
Visual attributes – Minor Defects	G-1	2.5
Physical Dimension	S-2 / N=13	4.0
Physical Properties	S-2	4.0
Force at break	N=13	N/A
Extractable Protein Content	N=3	N/A
Powder Amount	N=2	N/A

PACKAGING

Material	Quantity	Dimension (cm)
Inner Box	100 Pcs.	22(L) x 11(W) x 6.5(H)
Master Carton	1000 Pcs.	45.5(L) x 34.5(W) x 22.5(H)



Powdered Non Sterile Latex Examination Gloves