

# Surgicare®

## DHD UNDERGLOVE PREMIUM

### Sterile Latex Indicator UnderGlove Powderfree Surgical Gloves with Moisturizer

www.kanamlatex.com

Surgicare® Underglove Premium Latex Powderfree Surgical gloves with **moisturizer** offers perfect barrier protection whenever double gloving is required. Its dexterity and sensitivity proves to make it an ideal **Underglove**. These gloves are non-chlorinated and Antitack Polymer Treated Micro Rough Textured.

#### SPECIAL SAFETY FEATURES

- ▶ **Dark Green coloured primary glove for Double gloving and Puncture Detection of secondary glove.**
- ▶ **Damp Hand Donning (DHD) with Moisturizer.** DHD Technology assists in donning gloves over damp hands and moisturizer keeps the skin hydrated & soft
- ▶ Powderfree to eliminate risks associated with glove powder.
- ▶ Extra Sensitive, Microtextured & Anatomical.
- ▶ Smooth inner Surface with Polymer coating enables smooth and easy donning of glove.
- ▶ Low Protein content, no more than 50µg of protein per gram and dm<sup>2</sup>, provides safety to doctor and patient.
- ▶ Non-detectable levels of Chemical Residue, minimize risk of skin irritation.
- ▶ The thin gloves enables easy donning of secondary glove, providing added protection in High Risk Surgery.
- ▶ Can be used as an **Underglove** or Individually.
- ▶ Manufactured in dust free environment with low bioburden and therefore reduced pyrogen.
- ▶ Sterilized in a validated process cycle as per ISO 11135 norms SAL 10<sup>-6</sup>.

#### STANDARDS

- ▶ ASTM D 3577 & EN 455-1, EN 455-2, EN 455-3 & ISO 10282.
- ▶ CE Mark - MDR: Regulation (EU) 2017/745, Class IIa.

#### QUALITY CONTROL

- ▶ 100% Inspected.
- ▶ Factory Standard of AQL 0.65 upon request.
- ▶ Passes Biocompatibility test as per ISO 10993-10.
- ▶ Passes Viral Penetration Test - ASTM F 1671.
- ▶ Tested for Chemotherapy Drugs as per ASTM D 6978.

Recommended for:

**All Surgeries, especially those of High Risk**  
**Excellent as an UNDERGLOVE**



**Kanam**  
**KANAM LATEX**  
**INDUSTRIES PVT. LTD.**

The trusted name the world over

Ooppotttil Buildings, K.K. Road  
Kottayam-686 001, Kerala, India  
Tel: 0481-2560108, 2563513  
Fax: 0481-2563614

E-mail: [customercare@kanamlatex.com](mailto:customercare@kanamlatex.com)  
[ho@kanamlatex.com](mailto:ho@kanamlatex.com)